

NHP2 Protein, Human, Recombinant (His)

General Information

Synonyms:	NOLA2;Nucleolar protein family A member 2;snoRNP protein NHP2;H/ACA ribonucleoprotein complex subunit 2;NHP2;NHP2P.
Protein Construction:	Met1-Leu153
Species:	Human
Expression Host:	E. coli
Accession:	Q9NX24
Molecular Weight:	22 KDa (reducing condition)
AA Sequence:	Met1-Leu153

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 100 mM NaCl, 1 mM DTT, pH 8.0.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

NHP2 belongs to the H/ACA snoRNPs gene family. snoRNPs are involved in various aspects of rRNA processing and modification and have been organized into 2 families: C/D and H/ACA. NHP2 forms a small ribonucleoprotein particle with GAR1 (NOLA1). NHP2 is involved in a variety of aspects of rRNA processing and modification. NHP2 localizes to the dense fibrillar component of the nucleolus and in nuclear Cajal bodies. NHP2 may also be required for correct processing or intranuclear trafficking of TERC, the RNA component of the telomerase reverse transcriptase (TERT) holoenzyme.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481