

## Syntaxin 8 Protein, Human, Recombinant

### General Information

Synonyms:	STX8;Syntaxin-8
Protein Construction:	Met1-Gly215
Species:	Human
Expression Host:	E. coli
Accession:	Q9UNK0
Molecular Weight:	28 KDa (reducing condition)
AA Sequence:	Met1-Gly215

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 20 mM PB, 150 mM NaCl, pH 7.4.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Syntaxin-8 is a single-pass type IV membrane protein which belongs to the syntaxin family. It contains one t-SNARE coil homology domain. STX8 is highly expressed in heart, also found in brain, kidney, liver, lung, placenta, skeletal muscle, spleen and pancreas. STX8 is involved in protein trafficking from early to late endosomes via vesicle fusion and exocytosis. It as a vesicle trafficking protein functions in the early secretory pathway, possibly mediating retrograde transport form cis-golgi membrane to the ER.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481