

CXCL15 Protein, Mouse, Recombinant (His)

General Information

Synonyms:	Lungkine;C-X-C motif chemokine 15;Small-inducible cytokine B15;Scyb15;Cxcl15
Protein Construction:	Gln26-Ala167
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	Q9WVL7
Molecular Weight:	20-22 KDa (reducing condition)
AA Sequence:	Gln26-Ala167

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Mouse C-X-C motif chemokine 15, also known as Cxcl15, is a secreted protein which is member of the ELR motif-containing CXC chemokines. It expressed at low levels in fetal lung, the expression restricted to the lung, produced by bronchoepithelial cells and is released into the airways. It plays an important role in lung-specific neutrophil trafficking during normal and inflammatory conditions. It also appears chemotactic for neutrophils.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481