

NTNG1 Protein, Human, Recombinant (His)

General Information

Synonyms:	KIAA0976;Netrin-G1;Laminet-1;NTNG1;LMNT1
Protein Construction:	His29-Ser409
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q9Y2I2
Molecular Weight:	51 KDa (reducing condition)
AA Sequence:	His29-Ser409

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 20 mM PB, 150 mM NaCl, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Netrin-G1 (NTNG1) is a member of a conserved family of proteins that act as axon guidance cues during vertebrate nervous system development. Netrin-G1 contains one laminin EGF-like domain and one laminin N-terminal domain, Netrin-G1 is highly expressed in the thalamus, lowly in other tissue. Netrin-G1 localizes to the cell membrane. Netrin-G1 interacts with NGL1 and is glycosylated in the N-terminal. In addition, Netrin-G1 can promotes neurite outgrowth of both axons and dendrites.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481