

## LRP-12 Protein, Human, Recombinant (His)

### General Information

Synonyms:	LRP-12;ST7;Suppressor of Tumorigenicity 7 Protein;LRP12;Low-Density Lipoprotein Receptor-Related Protein 12
Protein Construction:	Glu33-Ile488
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q9Y561
Molecular Weight:	68 KDa (reducing condition)
AA Sequence:	Glu33-Ile488

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 20 mM PB, 150 mM NaCl, pH 7.2.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Low-Density Lipoprotein Receptor-Related Protein 12 (LRP12) belongs to the LDLR family. LRP12 is a type I transmembrane protein and widely expressed in heart, skeletal muscle, brain, lung, placenta and pancreas. LRP12 contains 2 CUB domain and 5 LDL-receptor class A domain. LRP12 has been shown to interact with GNB2L1, ZFYVE9 and ITGB1BP3. LRP12 is a receptor probably, which may be involved in the internalization of lipophilic

molecules and/or signal transduction. In addition, LRP12 may act as a tumor suppressor.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481