

BCMA/TNFRSF17 Protein, Cynomolgus, Recombinant (hFc)

General Information

Synonyms:	Bcma;Tnfrsf17;B-cell maturation protein;Bcm;CD269;Tumor necrosis factor receptor superfamily member 17
Protein Construction:	Met1-Ala53
Species:	Cynomolgus
Expression Host:	HEK293 Cells
Accession:	G7Q0I4
Molecular Weight:	38-45 KDa (reducing condition)
AA Sequence:	Met1-Ala53

QC Testing

Biological Activity:	1. Loaded Human APRIL-His on HIS1K Biosensor, can bind Cynomolgus BCMA-Fc with an affinity constant of 0.54 nM as determined in BLI assay. (Regularly tested) 2. Immobilized Human APRIL -Flag-His at 1µg/ml (100 µl/well) can bind Cynomolgus BCMA-Fc. The ED50 of Cynomolgus BCMA-Fc is 4.23 ng/ml. (Regularly tested)
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/µg (1 EU/µg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 50 mM Tris-HCl, 100 mM Glycine, pH 7.5.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

B cell maturation antigen (BCMA) is a member of the TNF receptor superfamily. It has been designated TNFRSF17. BCMA is a type III membrane protein containing one extracellular cysteine rich domain. Within the TNFRSF, it

A DRUG SCREENING EXPERT

shares the highest homology with TACI. BCMA and TACI have both been shown to bind to APRIL and BAFF, members of the TNF ligand superfamily. BCMA expression has been found in immune organs and mature B cell lines. Although some expression has been observed at the cell surface, BCMA appears to be localized to the Golgi compartment. The binding of BCMA to APRIL or BAFF has been shown to stimulate IgM production in peripheral blood B cells and increase the survival of cultured B cells. This data suggests that BCMA may play an important role in B cell development, function and regulation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481