

CD28 Protein, Human/Cynomolgus/Rhesus macaque, Recombinant (His & Avi), Biotinylated

General Information

Synonyms:	CD28 antigen (Tp44);CD28 molecule;CD28;MGC138290
Protein Construction:	Asn19-Pro152
Species:	Human,Cynomolgus,Rhesus
Expression Host:	HEK293 Cells
Accession:	P10747
Molecular Weight:	18.2 kDa (predicted). Due to glycosylation, the protein migrates to 38-50 kDa based on Tris-Bis PAGE result.

QC Testing

Biological Activity:	Immobilized Anti-CD28 Antibody, hFc Tag at 2 μ g/ml (100 μ l/well) on the plate. Dose response curve for Biotinylated Human/Cynomolgus/Rhesus macaque CD28, His Tag with the EC50 of 24.4ng/ml determined by ELISA.
Purity:	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS (pH 7.4). Typically, 8% trehalose is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

CD28 is the receptor for CD80 (B7-1) and CD86 (B7-2) proteins. When activated by Toll-like receptor ligands, the CD80 expression is upregulated in antigen presenting cells (APCs). The CD86 expression on antigen presenting cells is constitutive. CD28 is the only B7 receptor constitutively expressed on naive T cells.

Reference

Suntharalingam G, et al. Cytokine Storm in a Phase 1 Trial of the Anti-CD28 Monoclonal Antibody TGN1412[J]. New England Journal of Medicine, 2006, 355(10):1018-1028.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481