

PVRIG Protein, Mouse, Recombinant (hFc & Avi)

General Information

Synonyms:	PREL-2;RMO1;MGC2463;C7orf15;MGC104322;PREL2;PVRIG;LPD;MGC138297;RaIGDS/AF-6;ALS2CR9;CD112R;ALS2CR18
Protein Construction:	Ser35-Asp165
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	A0A1B0GS01
Molecular Weight:	42.70 kDa (predicted). Due to glycosylation, the protein migrates to 53-63 kDa based on Tris-Bis PAGE result.

QC Testing

Biological Activity:	Immobilized Mouse Nectin-2, His Tag at 2µg/ml (100µl/well) on the plate. Dose response curve for Mouse PVRIG, hFc Tag with the EC50 of 0.81ug/ml determined by ELISA.
Purity:	> 95% as determined by Tris-Bis PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS (pH 7.4). Typically, 8% trehalose is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Murine PVRIG interacted weakly with poliovirus receptor (PVR) but bound poliovirus receptor-like 2 (PVRL2) strongly, making the latter its principal ligand. PVRIG is an inducible checkpoint receptor and that targeting PVRIG-PVRL2 interactions results in increased CD8+ T-cell function and reduced tumor growth.

Reference

Murter B, et al. Mouse PVRIG Has CD8+ T Cell-Specific Coinhibitory Functions and Dampens Antitumor Immunity. Cancer Immunol Res. 2019 Feb;7(2):244-256. doi: 10.1158/2326-6066.CIR-18-0460. Epub 2019 Jan 18. PMID: 30659055; PMCID: PMC7038785.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481