

SIRP alpha V2 Protein, Human, Recombinant (aa 31-369, His)

General Information

Synonyms:	MYD-1;SHPS1;MFR;PTPNS1;CD172a;SIRPA;SIRP α V2/CD172a;SHPS-1;BIT;MYD1;SIRP α 2; SIRPalpha2;P84;SIRP α ;SIRP alpha
Protein Construction:	Glu31-Arg369
Species:	Human
Expression Host:	HEK293 Cells
Accession:	AAH26692.1
Molecular Weight:	38 kDa (predicted). Due to glycosylation, the protein migrates to 55-65 kDa based on Tris-Bis PAGE result.

QC Testing

Biological Activity:	1. Immobilized Human SIRP alpha V2, His Tag at 1 μ g/ml (100 μ l/Well) on the plate. Dose response curve for Human CD47, hFc Tag with the EC50 of 0.15 μ g/ml determined by ELISA. 2. Human SIRP alpha V2, His Tag captured on Protein A Chip via Anti-His Antibody can bind Anti-SIRP alpha Ab., hFc Tag with an affinity constant of 9.11nM as determined in SPR assay.
Purity:	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS (pH 7.4). Typically, 8% trehalose is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Signal regulatory protein α (SIRP α) is a regulatory membrane glycoprotein from SIRP family expressed mainly by myeloid cells and also by stem cells or neurons. SIRP α acts as inhibitory receptor and interacts with a broadly

A DRUG SCREENING EXPERT

expressed transmembrane protein CD47 also called the "don't eat me" signal. Cancer cells highly expressed CD47 that activate SIRP α and inhibit macrophage-mediated destruction.

Reference

Weiskopf K, et al. Engineered SIRP α variants as immunotherapeutic adjuvants to anticancer antibodies[J]. Science, 2013, 341(6141):88-91.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481