

## CD9P1 Protein, Mouse, Recombinant (His)

### General Information

Synonyms:	CD9P-1;EWI-F;CD9 partner 1;CD315
Protein Construction:	Arg22-Pro832
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	Q9WV91
Molecular Weight:	92.21 kDa (predicted). Due to glycosylation, the protein migrates to 110-120 kDa based on Tris-Bis PAGE result.

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 $\mu$ m filter, containing PBS (pH 7.4). Typically, 8% trehalose is incorporated as a protective agent before lyophilization.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

The membrane protein CD9P-1 is a major component of the tetraspanin web, a network of molecular interactions in the plasma membrane, in which it specifically associates with tetraspanins CD9 and CD81. All CD9P-1 isoforms associate with CD9 leading to additional level of complexity of this primary complex in the tetraspanin web.

Reference

André M, et al. Glycosylation status of the membrane protein CD9P- Proteomics. 2007 Nov;7(21):3880-95. doi: 10.1002/pmic.200700355. PMID: 17960739.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481