


KU-0058684

Chemical Properties

CAS No. : 623578-11-0
 Formula: C₁₉H₁₄FN₃O₃
 Molecular Weight: 351.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	KU-0058684 is a potent PARP and DNA-PK inhibitors.
Targets(IC50)	Others,PARP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8463 mL	14.2316 mL	28.4633 mL
5 mM	0.5693 mL	2.8463 mL	5.6927 mL
10 mM	0.2846 mL	1.4232 mL	2.8463 mL
50 mM	0.0569 mL	0.2846 mL	0.5693 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481