


## Dihydroisoferulic acid

## Chemical Properties

CAS No. :	1135-15-5
Formula:	C <sub>10</sub> H <sub>12</sub> O <sub>4</sub>
Molecular Weight:	196.20
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dihydroisoferulic acid is a useful organic compound for research related to life sciences. The catalog number is T126359 and the CAS number is 1135-15-5.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0968 mL	25.4842 mL	50.9684 mL
5 mM	1.0194 mL	5.0968 mL	10.1937 mL
10 mM	0.5097 mL	2.5484 mL	5.0968 mL
50 mM	0.1019 mL	0.5097 mL	1.0194 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481