

Surangin B

Chemical Properties


CAS No. :

Formula: C₂₉H₃₈O₇

Molecular Weight: 498.62

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Surangin B is a natural product and has a wide range of applications in life science related research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0055 mL	10.0277 mL	20.0554 mL
5 mM	0.4011 mL	2.0055 mL	4.0111 mL
10 mM	0.2006 mL	1.0028 mL	2.0055 mL
50 mM	0.0401 mL	0.2006 mL	0.4011 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481