

LEI-301

## Chemical Properties


CAS No. : 2571585-09-4

Formula: C<sub>25</sub>H<sub>25</sub>NO<sub>3</sub>

Molecular Weight: 387.47

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	LEI-301 is a novel potent inhibitor for the PLAAT family members, reducing the NAE levels, including anandamide, in cells overexpressing PLAAT2 or PLAAT5.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5808 mL	12.9042 mL	25.8084 mL
5 mM	0.5162 mL	2.5808 mL	5.1617 mL
10 mM	0.2581 mL	1.2904 mL	2.5808 mL
50 mM	0.0516 mL	0.2581 mL	0.5162 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481