


WAY-639889

Chemical Properties

CAS No. : 432506-24-6
 Formula: C₂₅H₂₇N₅O
 Molecular Weight: 413.51
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-639889 is a small molecule compound with selective and potent inhibitory effects on neuropeptide Y-5 receptors.
Targets(IC50)	Neuropeptide Y Receptor

Solubility Information

Solubility	DMSO: 50 mg/mL (120.92 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4183 mL	12.0916 mL	24.1832 mL
5 mM	0.4837 mL	2.4183 mL	4.8366 mL
10 mM	0.2418 mL	1.2092 mL	2.4183 mL
50 mM	0.0484 mL	0.2418 mL	0.4837 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Timothy J. et al. Hepatitis C Virus Entry Inhibitors. WO2008021745A2.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481