


2'-O-Methyl-2-thiouridine

Chemical Properties

CAS No. :	113886-72-9
Formula:	C10H14N2O5S
Molecular Weight:	274.29
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-O-Methyl-2-thiouridine is a purine nucleoside analogue present in synthetic thermophilic bacterial tRNAs. 2'-O-Methyl-2-thiouridine is more selective for A than unmodified U.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6458 mL	18.2289 mL	36.4578 mL
5 mM	0.7292 mL	3.6458 mL	7.2916 mL
10 mM	0.3646 mL	1.8229 mL	3.6458 mL
50 mM	0.0729 mL	0.3646 mL	0.7292 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Robak T, et al. Purine nucleoside analogs in the treatment of rarer chronic lymphoid leukemias. *Curr Pharm Des.* 2012;18(23):3373-88.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481