

## Benzyl-PEG2-ethanol

## Chemical Properties

CAS No. : 2055300-50-8


Formula: C<sub>13</sub>H<sub>20</sub>O<sub>3</sub>

Molecular Weight: 224.3

Storage: Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Benzyl-PEG2-ethanol, a PEG-based linker for PROTACs, joins two essential ligands crucial for forming PROTAC molecules and enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are composed of two distinct ligands linked together: one ligand targets an E3 ubiquitin ligase, while the other targets a specific protein. By harnessing the intracellular ubiquitin-proteasome system, PROTACs selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4583 mL	22.2916 mL	44.5831 mL
5 mM	0.8917 mL	4.4583 mL	8.9166 mL
10 mM	0.4458 mL	2.2292 mL	4.4583 mL
50 mM	0.0892 mL	0.4458 mL	0.8917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481