

RS 116 0086

## Chemical Properties


CAS No. : 148702-61-8

Formula: C<sub>19</sub>H<sub>27</sub>N<sub>2</sub>O<sub>4</sub>

Molecular Weight: 333.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	RS 116 0086 is a bioactive chemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9992 mL	14.9961 mL	29.9922 mL
5 mM	0.5998 mL	2.9992 mL	5.9984 mL
10 mM	0.2999 mL	1.4996 mL	2.9992 mL
50 mM	0.060 mL	0.2999 mL	0.5998 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481