

SR 26831

## Chemical Properties


CAS No. : 139995-67-8

Formula: C<sub>19</sub>H<sub>23</sub>ClNO<sub>3</sub>S

Molecular Weight: 380.91

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	SR 26831 is a selective elastase inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6253 mL	13.1265 mL	26.2529 mL
5 mM	0.5251 mL	2.6253 mL	5.2506 mL
10 mM	0.2625 mL	1.3126 mL	2.6253 mL
50 mM	0.0525 mL	0.2625 mL	0.5251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481