

Neocryptomerin

Chemical Properties

CAS No. :	20931-36-6
Formula:	C ₃₁ H ₂₀ O ₁₀
Molecular Weight:	552.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Neocryptomerin is a useful organic compound for research related to life sciences and the catalog number is T125374.
Targets(IC50)	Others,SARS-CoV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.810 mL	9.0499 mL	18.0999 mL
5 mM	0.362 mL	1.810 mL	3.620 mL
10 mM	0.181 mL	0.905 mL	1.810 mL
50 mM	0.0362 mL	0.181 mL	0.362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481