


Limocitrin 3-O-sophoroside

Chemical Properties

CAS No. :	114882-19-8
Formula:	C ₂₉ H ₃₄ O ₁₈
Molecular Weight:	670.57
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Limocitrin 3-O-sophoroside is a compound that can be used in biological research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4913 mL	7.4563 mL	14.9127 mL
5 mM	0.2983 mL	1.4913 mL	2.9825 mL
10 mM	0.1491 mL	0.7456 mL	1.4913 mL
50 mM	0.0298 mL	0.1491 mL	0.2983 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481