


Ispronicline hemigalactarate

Chemical Properties

CAS No. :	252870-54-5
Formula:	C ₁₄ H ₂₂ N ₂ O ₁₁ /2C ₆ H ₁₀ O ₈
Molecular Weight:	339.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ispronicline, also known as AZD3480; TC 1734; TC-0173, and RJR 1734, is a $\alpha_4\beta_2$ nicotinic acetylcholine receptor-selective agonist with cognition-enhancing and neuroprotective properties. Ispronicline is potentially useful for the treatment of attention-deficit/hyperactivity disorder in adults.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9463 mL	14.7314 mL	29.4629 mL
5 mM	0.5893 mL	2.9463 mL	5.8926 mL
10 mM	0.2946 mL	1.4731 mL	2.9463 mL
50 mM	0.0589 mL	0.2946 mL	0.5893 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481