

Fmoc-Gly3-Val-Cit-PAB

Chemical Properties


CAS No. : 2647914-09-6

Formula: C₃₉H₄₈N₈O₉

Molecular Weight: 772.85

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fmoc-Gly3-Val-Cit-PAB is a cleavable linker used in the assembly of antibody-drug conjugates (ADCs).
Targets(IC50)	ADC Linker
In vitro	ADCs (Antibody-Drug Conjugates) consist of an antibody that is connected to an ADC cytotoxin via an ADC linker [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2939 mL	6.4696 mL	12.9391 mL
5 mM	0.2588 mL	1.2939 mL	2.5878 mL
10 mM	0.1294 mL	0.647 mL	1.2939 mL
50 mM	0.0259 mL	0.1294 mL	0.2588 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481