

NHS-SS-Ph

## Chemical Properties


CAS No. : 2750799-12-1

Formula: C<sub>16</sub>H<sub>17</sub>N<sub>1</sub>O<sub>7</sub>S<sub>2</sub>

Molecular Weight: 399.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	NHS-SS-Ph, a cleavable linker utilized in antibody-drug conjugate (ADC) synthesis, facilitates the release of therapeutic agents under specific intracellular conditions [1].
Targets(IC50)	ADC Linker
In vitro	ADCs (antibody-drug conjugates) consist of an antibody connected to an ADC cytotoxin via an ADC linker.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5035 mL	12.5175 mL	25.035 mL
5 mM	0.5007 mL	2.5035 mL	5.007 mL
10 mM	0.2504 mL	1.2518 mL	2.5035 mL
50 mM	0.0501 mL	0.2504 mL	0.5007 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481