


Propargyl-PEG4-GGFG-DXd

Chemical Properties

CAS No. :	2762518-94-3
Formula:	C54H63FN8O15
Molecular Weight:	1083.12
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Propargyl-PEG4-GGFG-DXd is a drug-linker conjugate for antibody-drug conjugates (ADCs), incorporating an ADC linker and the DNA topoisomerase I inhibitor DXd [1].
Targets(IC50)	Drug-Linker Conjugates for ADC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9233 mL	4.6163 mL	9.2326 mL
5 mM	0.1847 mL	0.9233 mL	1.8465 mL
10 mM	0.0923 mL	0.4616 mL	0.9233 mL
50 mM	0.0185 mL	0.0923 mL	0.1847 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481