

SB1-G-187

## Chemical Properties

CAS No. : 2769753-18-4


Formula: C52H54F3N9O10

Molecular Weight: 1022.03

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	SB1-G-187 is a PROTAC (PROteolysis TArgeting Chimera) that functions as a multi-kinase degrader [1].
Targets(IC50)	PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9784 mL	4.8922 mL	9.7844 mL
5 mM	0.1957 mL	0.9784 mL	1.9569 mL
10 mM	0.0978 mL	0.4892 mL	0.9784 mL
50 mM	0.0196 mL	0.0978 mL	0.1957 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481