

## PROTAC IRAK3 degrade-1 formic

## Chemical Properties

CAS No. :

Formula: C48H65N11O8

Molecular Weight: 924.1

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PROTAC IRAK3 Degrade-1 (Compound 23) is a selective and potent IRAK3 degrader, demonstrating an IC50 value of 5 nM [1].
Targets(IC50)	PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0821 mL	5.4107 mL	10.8213 mL
5 mM	0.2164 mL	1.0821 mL	2.1643 mL
10 mM	0.1082 mL	0.5411 mL	1.0821 mL
50 mM	0.0216 mL	0.1082 mL	0.2164 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481