


Laminin (925-933)TFA

Chemical Properties

CAS No. :	2828432-44-4
Formula:	C42H63F3N12O16S
Molecular Weight:	1081.08
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Laminin (925-933) (TFA) is a peptide corresponding to the 925-933 sequence of the Laminin B1 chain, known to bind the laminin receptor [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.925 mL	4.625 mL	9.250 mL
5 mM	0.185 mL	0.925 mL	1.850 mL
10 mM	0.0925 mL	0.4625 mL	0.925 mL
50 mM	0.0185 mL	0.0925 mL	0.185 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481