

Arg-Gly-Glu-Ser TFA

Chemical Properties

CAS No. :

Formula: C18H30F3N7O10

Molecular Weight: 561.47

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Arg-Gly-Glu-Ser TFA, an RGD-related peptide, serves as a control for evaluating the inhibitory activity of RGDS on fibrinogen binding to activated platelets.
Targets(IC50)	Integrin
In vivo	The co-administration of Arg-Gly-Glu-Ser (5 mg/kg) with either LPS or saline, along with RGDS, does not influence the counts of neutrophils and macrophages, nor does it have an impact on protein accumulation, when compared to mice treated with LPS or saline alone [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.781 mL	8.9052 mL	17.8104 mL
5 mM	0.3562 mL	1.781 mL	3.5621 mL
10 mM	0.1781 mL	0.8905 mL	1.781 mL
50 mM	0.0356 mL	0.1781 mL	0.3562 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481