

WH-10417-099

## Chemical Properties

CAS No. : 2769753-42-4


Formula: C52H62N8O13S

Molecular Weight: 1039.16

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

|               |  |
|---------------|--|
| Description   | WH-10417-099 is a degrader capable of inducing the degradation of over 125 unique kinases [1]. |
| Targets(IC50) | Others,PROTACs   |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg      |
|-------|-----------|-----------|-----------|
| 1 mM  | 0.9623 mL | 4.8116 mL | 9.6232 mL |
| 5 mM  | 0.1925 mL | 0.9623 mL | 1.9246 mL |
| 10 mM | 0.0962 mL | 0.4812 mL | 0.9623 mL |
| 50 mM | 0.0192 mL | 0.0962 mL | 0.1925 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481