

## B-Raf IN 15

## Chemical Properties

CAS No. : 832107-31-0


Formula: C<sub>19</sub>H<sub>15</sub>N<sub>3</sub>O<sub>3</sub>

Molecular Weight: 333.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	B-Raf IN 15 is a BRAF inhibitor that inhibits BRAFWT and BRAFV600E and can be used in the study of melanoma and cancer, and can be optimized for use as a more effective BRA F inhibitor.
Targets(IC50)	Raf

## Solubility Information

Solubility	DMSO: 80 mg/mL (239.94 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 3.3 mg/mL (9.9 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9993 mL	14.9966 mL	29.9931 mL
5 mM	0.5999 mL	2.9993 mL	5.9986 mL
10 mM	0.2999 mL	1.4997 mL	2.9993 mL
50 mM	0.060 mL	0.2999 mL	0.5999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Qin J, et al. Identification of a novel family of BRAF(V600E) inhibitors. J Med Chem. 2012 Jun 14;55(11):5220-30.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481