


(R)-Azasetron besylate

Chemical Properties

CAS No. :	2025360-91-0
Formula:	C ₂₃ H ₂₆ ClN ₃ O ₆ S
Molecular Weight:	507.99
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(R)-Azasetron besylate (SENS-401), an orally active calcineurin inhibitor, has been demonstrated to mitigate cisplatin-induced hearing loss and cochlear damage [1][2].
Targets(IC50)	Phosphatase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9685 mL	9.8427 mL	19.6854 mL
5 mM	0.3937 mL	1.9685 mL	3.9371 mL
10 mM	0.1969 mL	0.9843 mL	1.9685 mL
50 mM	0.0394 mL	0.1969 mL	0.3937 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481