

## Anti-Mouse CD209b Antibody (22D1)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at -20°C  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse CD209b Antibody is an IgG1 inhibitor targeting mouse CD209b, sourced from the Armenian hamster.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481