

## Anti-Mouse CD3 Antibody (17A2)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse CD3 Antibody (17A2) is a murine IgG2b class antibody targeting CD3, which can be used to study TCR signaling, T cell activation and antigen recognition.
Targets(IC50)	Others,Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481