

## Anti-Mouse Ly-6G/Ly-6C Antibody (RB6-8C5)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse Ly-6G/Ly-6C Antibody is an IgG2b isotype antibody targeting Ly-6G and Ly-6C antigens on mouse cells, functioning as an inhibitor.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481