

## Alefacept

## Chemical Properties

CAS No. : 222535-22-0  
Formula:  
Molecular Weight:  
Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Alefacept (BG 9273), a fusion protein combining human lymphocyte function-associated antigen 3 with immunoglobulin 1, is utilized in the research of chronic plaque psoriasis [1].
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481