

Fazpilodemab

Chemical Properties

CAS No. : 2517935-02-1

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fazpilodemab (BFKB8488A) is a humanized agonistic bispecific antibody that targets fibroblast growth factor receptor 1c (FGFR1c) and Klotho β . It is utilized in research pertaining to Type 2 diabetes mellitus (T2DM) and nonalcoholic fatty liver disease (NAFLD) [1].
Targets(IC50)	FGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481