

Rocatinlimab

Chemical Properties

CAS No. : 2431972-52-8

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rocatinlimab represents a highly specialized humanized IgG1 monoclonal antibody specifically engineered to target the OX40 receptor (CD134) with exceptional binding affinity. It effectively facilitates the measurable blockade of the OX40-OX40L co-stimulatory axis while simultaneously inducing the targeted depletion of OX40-positive activated T cells through antibody-dependent cellular cytotoxicity (ADCC) across various preclinical experimental models to evaluate the restoration of immune homeostasis during strictly monitored laboratory observation periods.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481