

## Mouse IgG1 kappa, Isotype Control

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Mouse IgG1 kappa, Isotype Control is a mouse monoclonal antibody serving as an isotype control for mouse IgG1 $\kappa$ antibodies. Mouse IgG1 kappa, Isotype Control serves as a negative control to ensure more reliable assessment of antibody-specific binding within experimental systems. It also functions as a baseline reference in cell analysis, flow cytometry, immunofluorescence, immunoprecipitation, and other immunological assays, thereby enhancing the scientific rigor and reproducibility of data interpretation.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481