

## Mouse IgG2a kappa, Isotype Control

## Chemical Properties

CAS No. :

Formula:

Molecular Weight: 143.58 kDa

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Mouse IgG2a kappa, Isotype Control is an isotype antibody to murine-derived IgG2a kappa antibody.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.007 mL	0.0348 mL	0.0696 mL
5 mM	0.0014 mL	0.007 mL	0.0139 mL
10 mM	0.0007 mL	0.0035 mL	0.007 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481