

Human IgG1 (N297A) kappa, Isotype Control

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human IgG1 (N297A) kappa Isotype Control is a human-mouse chimeric monoclonal antibody serving as an isotype control specifically for the human IgG1 (N297A) κ antibody.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481