

Myelin Basic Protein(87-99) TFA

Chemical Properties

CAS No. :

Formula: C76H115N20O19F3

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Myelin Basic Protein (87-99) TFA, an encephalitogenic peptide, promotes basic protein-specific T cell proliferation and induces Th1 polarization in peripheral blood mononuclear cells, contributing to the pathogenesis of multiple sclerosis (MS) [1] [2].
Targets(IC50)	PKC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481