


Anti-amyloid agent-1

Chemical Properties

CAS No. :	1396778-73-6
Formula:	C ₂₁ H ₁₇ F ₃ N ₂ O ₃
Molecular Weight:	402.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Anti-amyloid agent-1 is a potent compound that inhibits amyloid aggregation, offering a promising approach for research into the treatment of amyloidosis [1].
Targets(IC50)	Beta Amyloid,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4853 mL	12.4264 mL	24.8527 mL
5 mM	0.4971 mL	2.4853 mL	4.9705 mL
10 mM	0.2485 mL	1.2426 mL	2.4853 mL
50 mM	0.0497 mL	0.2485 mL	0.4971 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481