

FLLRN

Chemical Properties

CAS No. : 141136-84-7


Formula: C₃₁H₅₁N₉O₇

Molecular Weight: 661.79

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	FLLRN is a biologically active peptide serving as a PAR1-specific antagonist.
Targets(IC50)	Protease-activated Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5111 mL	7.5553 mL	15.1105 mL
5 mM	0.3022 mL	1.5111 mL	3.0221 mL
10 mM	0.1511 mL	0.7555 mL	1.5111 mL
50 mM	0.0302 mL	0.1511 mL	0.3022 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481