

Carnegine (Pectenine, N-Methylsalsolidine)

Chemical Properties


CAS No. :

Formula: C₁₃H₁₉NO₂

Molecular Weight: 221.3

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | Carnegine (Pectenine, N-Methylsalsolidine) is a useful organic compound for research related to life sciences and the catalog number is T131396. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 4.5188 mL | 22.5938 mL | 45.1875 mL |
| 5 mM | 0.9038 mL | 4.5188 mL | 9.0375 mL |
| 10 mM | 0.4519 mL | 2.2594 mL | 4.5188 mL |
| 50 mM | 0.0904 mL | 0.4519 mL | 0.9038 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481