


## Chaetoglobosin C

## Chemical Properties

CAS No. : 50645-76-6  
 Formula: C<sub>32</sub>H<sub>36</sub>N<sub>2</sub>O<sub>5</sub>  
 Molecular Weight: 528.64  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Chaetoglobosin C (Compound 4), an anthraquinone-chromone derivative of Chaetomium globosum KMITL-N0802, exhibits anti-tuberculosis activity [1].
Targets(IC50)	Others,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8916 mL	9.4582 mL	18.9165 mL
5 mM	0.3783 mL	1.8916 mL	3.7833 mL
10 mM	0.1892 mL	0.9458 mL	1.8916 mL
50 mM	0.0378 mL	0.1892 mL	0.3783 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481