


COH1

## Chemical Properties

CAS No. : 20217-22-5  
 Formula: C<sub>11</sub>H<sub>10</sub>N<sub>2</sub>O<sub>3</sub>S  
 Molecular Weight: 250.27  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	COH1 is a ribonucleotide reductase (RR) inhibitor utilized in research pertaining to cancer, mitochondrial diseases, and neurodegenerative disorders [1].
Targets(IC50)	DNA/RNA Synthesis

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9957 mL	19.9784 mL	39.9568 mL
5 mM	0.7991 mL	3.9957 mL	7.9914 mL
10 mM	0.3996 mL	1.9978 mL	3.9957 mL
50 mM	0.0799 mL	0.3996 mL	0.7991 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481