

AChE/BChE-IN-15

Chemical Properties

CAS No. :

Formula: C₂₉H₃₀N₆O₃

Molecular Weight: 510.59

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AChE/BChE-IN-15 (Compound 6d) is an inhibitor of acetylcholinesterase (AChE) and butyrylcholinesterase (BChE), with IC ₅₀ values of 20 nM for AChE and 220 nM for BChE. It targets the catalytic anionic site (CAS) and the peripheral anionic site (PAS) within the enzymes' active centers. AChE/BChE-IN-15 is applicable in Alzheimer's disease research [1].
Targets(IC50)	Others,Cholinesterase (ChE)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9585 mL	9.7926 mL	19.5852 mL
5 mM	0.3917 mL	1.9585 mL	3.917 mL
10 mM	0.1959 mL	0.9793 mL	1.9585 mL
50 mM	0.0392 mL	0.1959 mL	0.3917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481