

HPK1-IN-39

## Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>27</sub>N<sub>7</sub>O<sub>2</sub>

Molecular Weight: 469.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HPK1-IN-39 (Compound 10n), a selective inhibitor of Hematopoietic Progenitor Kinase 1 (HPK1) with an inhibitory concentration (IC <sub>50</sub> ) of 29 nM, impedes the phosphorylation of SLP76. This compound is applicable in the research of cancer immunotherapy [1].
Targets(IC <sub>50</sub> )	Others,MAPK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1297 mL	10.6487 mL	21.2974 mL
5 mM	0.4259 mL	2.1297 mL	4.2595 mL
10 mM	0.213 mL	1.0649 mL	2.1297 mL
50 mM	0.0426 mL	0.213 mL	0.4259 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481