

MAO-B-IN-24

Chemical Properties

CAS No. :

Formula: C₂₀H₁₈N₂O

Molecular Weight: 302.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MAO-B-IN-24 (compound 11h) is a selective, reversible, competitive inhibitor of MAO-B with an IC ₅₀ value of 1.60 μM and an IC ₅₀ of 22.42 μM for MAO-A. At 10 μM concentration, it reduces AChE and BChE activities to 54.58% and 88.43%, respectively [1].
Targets(IC ₅₀)	Cholinesterase (ChE), Monoamine Oxidase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3072 mL	16.536 mL	33.0721 mL
5 mM	0.6614 mL	3.3072 mL	6.6144 mL
10 mM	0.3307 mL	1.6536 mL	3.3072 mL
50 mM	0.0661 mL	0.3307 mL	0.6614 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481